

The background of the slide is a traditional marbled paper pattern, featuring swirling veins of orange, teal, and cream. A light gray semi-circular frame is centered on the slide, containing the text.

BRAIN TUMOR CLASSIFICATION

and assorted issues.

GOAL

- Build a classifier that is resistant to noise under a variety of conditions.


CONDITIONS

- Types of noise
- Resolutions
 - Blur

EXPERIMENTAL DESIGN

1. Build a baseline classifier
2. Evaluate performance with "perfect" images
 3. Add condition to image
 4. Evaluate performance
5. Update classifier to optimize performance

Data: Kaggle, MRI images



Train Baseline Classifier

Add noise

**Evaluate
Performance on
Baseline
Classifier**

Vary Resolution

**Evaluate
Performance on
Baseline
Classifier**

Add Motion Blur

**Evaluate
Performance on
Baseline
Classifier**


Update Classifier to Optimize for Condition

DATA


- Data is from Kaggle
- Brain MRI Images showing 3 types of tumors and a “no tumor” category.
- Images from a variety of angles

	Train Quantity	Test Quantity
Pituitary	1457	300
Glioma	1321	300
Meningioma	1339	306
No Tumor	1595	405

EXAMPLE DATA



BASELINE CLASSIFIER EVALUATION



BASELINE CLASSIFIER EVALUATION


Class Labels

0: Glioma

1: Meningioma

2: No Tumor


3: Pituitary




WHAT IS GAUSSIAN NOISE?

- Gaussian noise refers to random variations in the pixel values of the image that follow a Gaussian distribution.
- These variations can be caused by factors such as electronic noise in the MRI machine or imperfections in the signal acquisition process.

ADDING GAUSSIAN NOISE.
STANDARD DEVIATION = 0.1
MEAN = 0




GAUSSIAN NOISE: BEFORE TRAINING



Accuracy = 31.5%

GAUSSIAN NOISE: AFTER RETRAINING



Accuracy = 93.83%

Class Labels

0: Glioma

1:

Meningioma

2: No Tumor


3: Pituitary

WHAT IS POISSON NOISE?


- Poisson noise arises from the statistical nature of photon or particle detection, and it can be especially pronounced in low light conditions.
- Poisson noise can result in a "salt-and-pepper" effect in images, where random bright or dark pixels are introduced.
- In MRI, Poisson noise can be caused by various factors, such as variations in magnetic field strength or in the number of spins that emit a signal and can be amplified using higher magnetic fields.

ADDING POISSON NOISE.

$LAM = 0.1$



POISSON NOISE: BEFORE TRAINING




Accuracy = 30.89%

Class Labels

0: Glioma
1: Meningioma
2: No Tumor
3: Pituitary

POISSON NOISE: AFTER RETRAINING




Accuracy = 94.87%


WHAT IS QUANTIZATION NOISE?

- Quantization noise refers to the error introduced when representing continuous analog signals in a discrete form.
- Quantization noise can occur during the digitization process where the continuous image signal is sampled and then converted to a digital representation with a limited number of bits.
- This can result in a loss of information and the introduction of noise in the form of rounding errors and signal approximation.

ADDING QUANTIZATION NOISE. 32 LEVELS NOISE LEVEL 10



QUANTIZATION NOISE: BEFORE TRAINING



Accuracy = 30.89%

Class Labels

0: Glioma


1:

Meningioma

2: No Tumor

3: Pituitary

QUANTIZATION NOISE: AFTER RETRAINING




Accuracy = 76.15%


WHAT IS MOTION BLUR?

- Motion blur may occur due to patient motion during image acquisition or scanner instability.
- This can result in image blurring and distortion, reducing the quality of the image and potentially affecting the accuracy of diagnostic information.

ADDING MOTION BLUR.



MOTION BLUR: BEFORE TRAINING




Accuracy = 30.89%

Class Labels

0: Glioma
1: Meningioma
2: No Tumor
3: Pituitary

MOTION BLUR: AFTER RETRAINING




Accuracy = 89.58%


WHAT IS LOWERING RESOLUTION?

- Lowering the resolution can be used to reduce the size of the image and make it easier to store and process.
- However, it can also result in a loss of important details and features that may be important for diagnosis or analysis.

LOWERING RESOLUTION BY A FACTOR OF 2



LOWER RESOLUTION: BEFORE TRAINING



Accuracy = 31.27%

Class Labels

0: Glioma


1:

Meningioma

2: No Tumor

3: Pituitary

LOWER RESOLUTION: AFTER RETRAINING



Accuracy = 77.68%

CONCLUSION

- Gaussian noise and Poisson noise did not affect the performance very much after the model was retrained.
- Quantization noise suffered, even after retraining (76%).
 - Probably due to the magnitude of noise added.
- Motion blur maintained almost a 90% accuracy.
 - Despite the high level of blur added.
- Lowering the resolution also lowered the accuracy (77%).
 - Probably due to the level at which the resolution was lowered.